

B-8505 Stain Resistant White Slide Label

TDS No. B-8505
Effective Date: 02/27/2012

Description:
B-8505 is a stain resistant white polyester slide label

GENERAL

Print Technology: Thermal transfer
Material Type: Polyester
Finish: Matte white
Adhesive: Permanent acrylic
Liner: 55 lb. silicone release glassine

APPLICATIONS

Slide labeling and other lab identification

RECOMMENDED RIBBONS

Brady Series R4300

SPECIAL FEATURES

B-8505, when printed with the R4300 thermal transfer ribbon, is designed for use in the slide staining process and other laboratory identification which requires superior chemical resistance.

Details:

PHYSICAL PROPERTIES	TEST METHODS	AVERAGE RESULTS
Thickness	ASTM D 1000 -Total (excluding liner)	0.0038 inches (0.098 mm)
Adhesion to:	20 minute dwell 24 hour dwell	50 oz/inch (55.2 N/100mm) 55 oz/inch (59.7 N/100mm)
-Glass	20 minute dwell 24 hour dwell	73 oz/inch (80.3 N/100mm) 81 oz/inch (88.3 N/100mm)
-Polypropylene	20 minute dwell 24 hour dwell	48 oz/inch (52.8 N/100mm) 80 oz/inch (87.7 N/100mm)
-Stainless steel	20 minute dwell 24 hour dwell ASTM D 2979	48 oz/inch (52.8 N/100mm) 80 oz/inch (87.7 N/100mm)
Tack	Polyken™ Probe Tack 1 second dwell	37oz (1038 grams)

PERFORMANCE PROPERTIES-ENVIRONMENTAL RESISTANCE

B-8505 samples were printed with the R4300 Series ribbon on Bradyprinter™ THT Model 300X-Plus thermal transfer printer. Printed samples were laminated to glass microscope slides and then exposed to the indicated environmental conditions.

PERFORMANCE PROPERTIES	TEST METHODS	EFFECT TO LABEL	EFFECT TO PRINT IMAGE
High Service Temperature	30 days at 266°F (130°C)	No visible effect	No visible effect
Low Service Temperature	30 days at -112°F (-80°C)	No visible effect	No visible effect
Humidity Resistance*	30 days at 100°F (37°C), 80% relative humidity	No visible effect	No visible effect
Sterilizer resistance (autoclave resistance)	Ramp time 28-30 minutes, peak temperature 125°C, peak pressure 21psi and peak duration 5 minutes, Sterilemax Table Top Steam Sterilizer, Series 1277	No visible effect	No visible effect
Abrasion Resistance	Taber Abraser, CS-10 grinding wheels, 500 g/arm, 150 cycles	No visible effect	Print still legible after 500 cycles

PERFORMANCE PROPERTIES-CHEMICAL RESISTANCE

B-8505 samples were printed with R4300 Series ribbon on Bradyprinter™ THT Model 300X-Plus thermal transfer printer. Printed samples were laminated to aluminum panels and allowed to dwell 24 hours prior to testing. Testing consisted of 5 cycles of 10 minute immersions in the specified chemical followed by 30 minute recovery periods. After the final immersion, the printed image was rubbed 10 times with a cotton swab saturated with the test fluid. Testing was conducted at room temperature.

CHEMICAL REAGENT	EFFECT TO PRINT/TOPCOAT WITHOUT RUB	EFFECT TO PRINT/TOPCOAT WITH RUB	EFFECT TO ADHESIVE
Acetone	1	2	1
Ethanol	1	1	1
Toluene	1	2	1
Isopropanol	1	1	1
Xylene	1	2	1
MEK	1	2	1
Deionized Water	1	1	1
Dimethylsulfoxide (DMSO)	1	1	1
3% Alconox® Detergent	1	1	1
50% Acetic Acid	1	1	1
10% Sulfuric Acid	1	1	1
10% Clorox® Bleach Solution	1	1	1
10% Sodium Hydroxide	1	1	1

Rating scale:

1=no visible effect

2=slight smear or print removal or adhesive oozing, detectable but minimal

3=moderate smear or print removal or adhesive oozing

4=severe smear or print removal or adhesive oozing

5=complete smear, print removal or topcoat removal (delamination) or adhesive oozing

PERFORMANCE PROPERTIES-Slide Staining Test

Samples of B-8505 printed with R4300 Series ribbon were laminated to glass microscope slides. These slides were placed in the Leica AutoStainer XL automated slide stainer. The slides were processed using the H&E (hematoxilin and eosin) staining process.

Labels were adhered to the following brands of slides:

Cole-Parmer® - precleaned frosted, precleaned Superfrost® and plain precleaned Superfrost®

Fisher Scientific – precleaned frosted beveled

Erie Scientific - precleaned Superfrost® plus

Mercedes Medical - (precleaned, beveled edge, double frost)

Gold Seal Products® - Gold Seal® precleaned

REAGENT	STATION	STEP	DURATION (minutes)
None	Oven	1	10:00
Xylene	1	2	2:00
Xylene	2	3	2:00
100% denatured ethanol	3	4	2:00
100% denatured ethanol	4	5	2:00
70% denatured ethanol	5	6	1:00
Wash (tap water)	Wash 1	7	2:00
Hematoxilin stain (Harris formula modified, supplied by Ricca Chemical Company)	6	8	5:00
Wash (tap water)	Wash 2	9	2:00
Acid Alcohol (mixture of 1.25% of concentrated HCl, 69.26% of ethanol and 29.49% of DI H2O)	7	10	0.02
Wash (tap water)	Wash 3	11	3:00
Scott's bluing reagent (supplied by Fisher Scientific)	8	12	3:00
Wash (tap water)	Wash 4	13	3:00
Eosin Y stain (Stock formula, 1% alcoholic solution supplied by Ricca Chemical Company)	14	14	2:00
95% denatured ethanol	15	15	0:30
100% denatured ethanol	16	16	2:00
100% denatured ethanol	17	17	2:00

100% denatured ethanol	18	18	2:00
Xylene	Exit	19	30:00

**Test results: All the labels remained adhered to all types of slides through the process. Print was not affected.
Label surface was slightly pink after the process.**

Product testing, customer feedback, and history of similar products, support a customer performance expectation of at least **two years from the date of receipt** for this product as long as this product is stored in its original packaging in an environment *below 80 degrees F (27°C) and 60% RH*. We are confident that our product will perform well beyond this time frame. However, it remains the responsibility of the user to assess the risk of using such product. We encourage customers to develop functional testing protocols that will qualify a product's fitness for use, in their actual applications.

Trademarks:

ASTM: American Society for Testing and Materials (U.S.A.)
 Alconox® is a registered trademark of Alconox Co.
 BradyPrinter™ is a trademark of Brady Worldwide, Inc.
 Clorox® is a registered trademark of The Clorox Company
 Cole-Parmer® is a registered trademark of the Cole-Parmer Instrument Company
 Gold Seal® is a registered trademark of Erie Scientific Company
 Polyken™ is a trademark of Testing Machines Inc.
 Superfrost® is a registered trademark of Erie Scientific Company

Note: All values shown are averages and should not be used for specification purposes.

Test data and test results contained in this document are for general information only and shall not be relied upon by Brady customers for designs and specifications, or be relied on as meeting specified performance criteria. Customers desiring to develop specifications or performance criteria for specific product applications should contact Brady for further information.

Product compliance information is based upon information provided by suppliers of the raw materials used by Brady to manufacture this product or based on results of testing using recognized analytical methods performed by a third party, independent laboratory. As such, Brady makes no independent representations or warranties, express or implied, and assumes no liability in connection with the use of this information.

WARRANTY

Brady products are sold with the understanding that the buyers will test them in actual use and determine for themselves their adaptability to their intended uses. Brady warrants to the buyers that its products are free from defects in material and workmanship, but limits its obligation under this warranty to replacement of the product shown to Brady's satisfaction to have been defective at the time Brady sold it. This warranty does not extend to any persons obtaining the product from the buyers. This warranty is in lieu of any other warranty, express or implied, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose, and of any other obligations or liability on Brady's part. Under no circumstances will Brady be liable for any loss, damage, expense, or consequential damages of any kind arising in connection with the use, or inability to use, Brady's products.

Copyright 2012 Brady Worldwide, Inc. | All Rights Reserved
 Material may not be reproduced or distributed in any form without written permission.

Brady North America | 6555 W. Good Hope Rd | Milwaukee, WI 53223 | USA | Tel: 414-358-6600 | Fax: 800-292-2289

Trademarks:

ASTM: American Society for Testing and Materials (U.S.A.)
 Alconox® is a registered trademark of Alconox Co.
 BradyPrinter™ is a trademark of Brady Worldwide, Inc.
 Clorox® is a registered trademark of The Clorox Company
 Cole-Parmer® is a registered trademark of the Cole-Parmer Instrument Company
 Gold Seal® is a registered trademark of Erie Scientific Company
 Polyken™ is a trademark of Testing Machines Inc.
 Superfrost® is a registered trademark of Erie Scientific Company

Note: All values shown are averages and should not be used for specification purposes.

Test data and test results contained in this document are for general information only and shall not be relied upon by Brady customers for designs and specifications, or be relied on as meeting specified performance criteria. Customers desiring to develop specifications or performance criteria for specific product applications should contact Brady for further information.

Product compliance information is based upon information provided by suppliers of the raw materials used by Brady to manufacture this product or based on results of testing using recognized analytical methods performed by a third party,

independent laboratory. As such, Brady makes no independent representations or warranties, express or implied, and assumes no liability in connection with the use of this information.

WARRANTY

Brady products are sold with the understanding that the buyers will test them in actual use and determine for themselves their adaptability to their intended uses. Brady warrants to the buyers that its products are free from defects in material and workmanship, but limits its obligation under this warranty to replacement of the product shown to Brady's satisfaction to have been defective at the time Brady sold it. This warranty does not extend to any persons obtaining the product from the buyers. This warranty is in lieu of any other warranty, express or implied, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose, and of any other obligations or liability on Brady's part. Under no circumstances will Brady be liable for any loss, damage, expense, or consequential damages of any kind arising in connection with the use, or inability to use, Brady's products.

Copyright 2016 Brady Corporation Asia Pte. Ltd. | All Rights Reserved
Material may not be reproduced or distributed in any form without written permission.

Brady Asia | 1 Kaki Bukit Crescent | Singapore 416236 | Singapore | Tel: 65 6477.7261 | Fax: 65 6748.7248